

## Skinergy Plus Closed / Split Toe


Manufactured using the latest innovation in silicone technology, Skinergy Plus provides an extremely high level of cosmesis for the below-knee user, incorporating a uniquely seamless finish.

Easy to fit, enhanced toughness and improved realism makes Skinergy Plus our most impressive leg cover yet.



### Features and Benefits

- Integrated reinforcing mesh and advanced silicone material for a more durable finish
- Now available in 19 colour shades with unique TrueFinish™ micropigmentation and seamless finish for authentic realism
- Split or closed toe options to suit any foot size 220mm to 300mm
- Two calf sizes
- Easily removable for hardware adjustment, maintenance or repair
- Suitable for ankle disarticulation and partial foot applications
- Easy to clean using soap and water
- Colour swatch is available to order to demonstrate all 19 shades
- 3-month warranty



Semi-custom nail finish



Closed toe option

### Ordering Information

Part Number	Part Number	Foot Size	Ankle Circ (cm)		Mid Calf Circ (cm)		Top of Cover Circ (cm)		Cover Length From Floor to: (cm)	
			Min	Max	Min	Max	Min	Max	Mid Calf	Cover Top
Split Toe	Closed Toe	Size	Min	Max	Min	Max	Min	Max	Mid Calf	Cover Top
TLC22LS/E*	TNS22LS/E*	22-23	20	25	28	35	26	32.5	32	55
TLC22RS/E*	TNS22RS/E*	22-23	20	25	28	35	26	32.5	32	55
TLC22LL/E*	TNS22LL/E*	22-23	20	25	32	40	30	37.5	32	55
TLC22RL/E*	TNS22RL/E*	22-23	20	25	32	40	30	37.5	32	55
TLC24LS/E*	TNS24LS/E*	24-25	22	27.5	31	39	29	36	34.5	57
TLC24RS/E*	TNS24RS/E*	24-25	22	27.5	31	39	29	36	34.5	57
TLC24LL/E*	TNS24LL/E*	24-25	22	27.5	35	44	33	41	34.5	57
TLC24RL/E*	TNS24RL/E*	24-25	22	27.5	35	44	33	41	34.5	57
TLC26LS/E*	TNS26LS/E*	26-27	22	28	32	40	30	37.5	37	61
TLC26RS/E*	TNS26RS/E*	26-27	22	28	32	40	30	37.5	37	61
TLC26LL/E*	TNS26LL/E*	26-27	22	28	37	46	35	44	37	61
TLC26RL/E*	TNS26RL/E*	26-27	22	28	37	46	35	44	37	61
TLC28LS/E*	TNS28LS/E*	28-30	24	30	34	42.5	32	40	39.5	66
TLC28RS/E*	TNS28RS/E*	28-30	24	30	34	42.5	32	40	39.5	66
TLC28LL/E*	TNS28LL/E*	28-30	24	30	39	49	37	46	39.5	66
TLC28RL/E*	TNS28RL/E*	28-30	24	30	39	49	37	46	39.5	66

E\* denotes colour. For example, TNS22LS/E5.

Warranty: 3 months

### Technical Information

Part Number	Description
NCT9740	Skinergy Plus True Finish 25mm Punchtool
E24655	Skinergy Plus True Finish Swatch Kit
E26625/E*	Skinergy Plus True Finish Patch Kit
E26626/E*	Skinergy Plus True Finish Spare Kit

Skinergy covers can be requested with the following non-standard finishes: uncut, unpainted toenails. If a non-standard finish or bilateral cover is required, please specify when ordering.

### TrueFinish Colours



Colour Matching:

These swatches are printed as guidance only. For colour matching, the swatch kits should be used (part number E24655).