



Quick Grip Pliers

Provided with serrated jaws for a firm grip on small items such as nails, screws, nuts, bolts. Two different shaped apertures enabling hand tools with different diameter handles to be used. These pliers may be locked manually, allowing the user to relax tension on the operating system. The moving jaw also has a shaped section for pushing or holding against a bench.

Part Number	Description	Left / Right	Thread	Length	Weight
1 B18220	Quick Grip Pliers - Adult	Left	3/8" BSW	140mm (5 1/2")	283g (10 oz)
B05218	Quick Grip Pliers - Adult	Right	3/8" BSW	140mm (5 1/2")	283g (10 oz)
B09185	Replacement Spring	Left or Right	n/a	n/a	n/a

Tweezers

Suitable for fine precise work. The jaws are serrated and curved, for ease of use and clear visibility.

Part Number	Description	Left / Right	Thread	Length	Weight
2 B01024	Tweezers	Left	3/8" BSW	95mm (3 3/4")	57g (2 oz)
B01025	Tweezers	Right	3/8" BSW	95mm (3 3/4")	57g (2 oz)
B09439	Replacement Rubber Ring for Tweezers	n/a	n/a	n/a	n/a

Long Nose Pliers

Provided with tapered jaws for precise work and two different sized apertures for gripping. Jaws may be locked with thumb screw.

Part Number	Description	Left / Right	Thread	Length	Weight
3 D00172	Long Nose Pliers - Adult	Left or Right	3/8" BSW	120mm (4 3/4")	255g (9 oz)

Universal Tool Holder

Designed for a stable grip on tools, bars, rods and tubes. Jaws may be manually locked on to an object allowing the user to relax tension on the operating system.

Part Number	Description	Left / Right	Thread	Length	Weight
4 D00175	Universal Tool Holder - Adult	Left or Right	3/8" BSW	120mm (4 3/4")	340g (12 oz)

Note: A full range of adapter plates are available. Please see page 111.
Thread locking compound must be used when assembling adapters to functional devices.