



Quick Disconnect Wrist Unit

Part Number	Description	Hand Size	Height	Weight
5 QDAWA-40	40mm Wrist Assembly. Includes Lamination Ring, Coupling Piece and Coaxial Plug	6 ¾	20mm (¾")	38g (1 ⅜ oz)
QDAWA-45	45mm Wrist Assembly. Includes Lamination Ring, Coupling Piece and Coaxial Plug	7 ¼	25mm (1")	38g (1 ⅜ oz)
QDAWA-50	50mm Wrist Assembly. Includes Lamination Ring, Coupling Piece and Coaxial Plug	7 ¾	25mm (1")	42g (1 ½ oz)
QDAWA-54	54mm Wrist Assembly. Includes Lamination Ring, Coupling Piece and Coaxial Plug	8 ¼	25mm (1")	44g (1 ½ oz)

Quick Disconnect Wrist - Spares

Part Number	Description	Hand Size	Height
6 QDALR-40	40mm Lamination Ring	6 ¾	20mm (¾")
QDALR-45	45mm Lamination Ring	7 ¼	23mm (⅞")
QDALR-50	50mm Lamination Ring	7 ¾	23mm (⅞")
QDALR-54	54mm Lamination Ring	8 ¼	23mm (⅞")
7 QDACP	Coupling Piece. Suitable for 40, 45, 50 and 54mm Quick Disconnect Wrist. Includes lock ring.		
8 QDACO	Coaxial Plug. Suitable for 40, 45, 50 and 54mm Quick Disconnect Wrist. Includes lock ring and screw.		

Battery

Part Number	Description	Chemistry	Capacity	Dimensions	Weight
9 BLPA72	7.2V Lithium Polymer Battery	Lithium Polymer	800mAh	68 x 33 x 17 mm (2 ⅝" x 1 ¼" x ⅝")	45g (1 ½ oz)

Battery Charger

Part Number	Description	Region
10 LPDC01	Lithium Polymer Battery Desktop Charger	World

Battery Receptacles

Forearm mounted battery receptacles are supplied in kits with all necessary formers

Part Number	Description	Battery
11 BHQDA72	Lithium Polymer Battery Holder with 2-Pin Output	BLPA72

Note: Suitable for connection to Quick Disconnect Wrist connection

12 BHA72	Lithium Polymer Battery Holder with Locking Connector	BLPA72
----------	---	--------

Note: Suitable for connection to Connector Block Cable CBBHA72.