



**Body-Powered Gripper**

Curved jaws allow spherical and cylindrical objects to be held firmly. Two levels of grip can be selected by simple user control. Hard tips for fine control. This item requires a hand connection plate, please see pages 111.

Part Number	Description	Left / Right	Base	Length	Weight
1 D11377	Body-Powered Gripper - Adult	Left	2 ½" Diameter	140mm (5 ½")	200g (7 oz)
D11378	Body-Powered Gripper - Adult	Right	2 ½" Diameter	140mm (5 ½")	200g (7 oz)
2 E17996	Replacement Jaw Lining	Left	n/a	n/a	n/a
E17997	Replacement Jaw Lining	Right	n/a	n/a	n/a

**Note:** A full range of adapter plates are available. Please see page 111. Thread locking compound must be used when assembling adapters to functional devices.