

**Body-Powered Gripper**

Curved jaws allow spherical and cylindrical objects to be held firmly. Two levels of grip can be selected by simple user control. Hard tips for fine control.

Hand connection plates available separately, please see page 137.



Part Number	Description	Left/Right	Base	Length	Weight	
1	D11377	Body-Powered Gripper - Adult	Left	2 1/8" Diameter	140mm (5 1/2")	200g (7 oz)
	D11378	Body-Powered Gripper - Adult	Right	2 1/8" Diameter	140mm (5 1/2")	200g (7 oz)
2	E17996	Replacement Jaw Lining	Left	N/A	N/A	N/A
	E17997	Replacement Jaw Lining	Right	N/A	N/A	N/A



**Universal Tool Holder**

Designed for a stable grip on tools, bars, rods and tubes. Jaws may be manually locked on to an object allowing the user to relax tension on the operating system.

Part Number	Description	Left/Right	Thread	Length	Weight
D00175	Universal Tool Holder - Adult	Left or Right	3/8" BSW	145mm (5 3/4")	347g (12 1/4 oz)



**Light Tool Holder**

An adjustable clamp with two jaws for gripping around handled tools. The head may be set at various angles and jaws tightened by means of a wing nut.

Part Number	Description	Left/Right	Thread	Length	Weight
B05191	Light Tool Holder - Adult	Left or Right	3/8" BSW	130mm (5 1/8")	250g (9 oz)