

**Technical Information: Child Cable Operated and Spring Hand**

Size	Width at Knuckles	Cable Pull	Max Opening	Weight	Base to Thumb Tip
5 1/2	50mm (2")	29mm (1 1/8")	44mm (1 3/4")	84g (3oz)	76mm (3")
6	57mm (2 1/4")	25mm (1")	50mm (2")	130g (4 1/2oz)	102mm (4")



**Child Cable Operated Hand Front Pull - Voluntary Opening**

Child's Cable Operated Hand, with a front pull and voluntary opening. These hands are constructed using rugged plastic mouldings to withstand the demands of a growing child, with functional benefits and an improved cosmetic appearance.

**Note:** Round fixed B05441 wrist (45mm) is not compatible with **size 6** hands. Recommended use is with friction housing (p.143) and hand plate (p.143).

Part Number	Size	UK Size	Left/Right	Application	Wrist
B08361	5 1/2	2" x 3"	Left	Child 2-5 years	40mm (1 1/2")
B08360	5 1/2	2" x 3"	Right	Child 2-5 years	40mm (1 1/2")
C07331	6	2 1/4" x 4"	Left	Child 5-8 years	45mm (1 3/4")
C07330	6	2 1/4" x 4"	Right	Child 5-8 years	45mm (1 3/4")



**Child Cable Operated Hand Back Pull - Voluntary Opening**

Child's Cable Operated Hand, with a back pull and voluntary opening. These hands are constructed using rugged plastic mouldings to withstand the demands of a growing child, with functional benefits and an improved cosmetic appearance.

**Note:** Round fixed B05441 wrist (45mm) is not compatible with **size 6** hands. Recommended use is with friction housing (p.143) and hand plate (p.143).

Part Number	Size	UK Size	Left/Right	Application	Wrist
B08363	5 1/2	2" x 3"	Left	Child 2-5 years	40mm (1 1/2")
B08362	5 1/2	2" x 3"	Right	Child 2-5 years	40mm (1 1/2")
C07931	6	2 1/4" x 4"	Left	Child 5-8 years	45mm (1 3/4")
C07932	6	2 1/4" x 4"	Right	Child 5-8 years	45mm (1 3/4")

## Cable Operated and Spring Hands - **Child Hand**

Upper Limb Prosthetics  
Ordering Information  
Myoelectric Hands  
Myoelectric Elbows  
Myoelectric Power & Inputs  
Cable Operated & Spring Hands  
Functional Devices  
Cosmetic Hands & Gloves  
Modular Arm Systems  
Wrist Units  
Mech. Elbow Units & Joints  
Suspension & Control Systems  
Materials & Tools



### Child Spring Hand

Child's manually operated Spring Hand with voluntary opening. These hands are constructed using rugged plastic mouldings to withstand the demands of a growing child, with functional benefits and an improved cosmetic appearance.

**Note:** Round fixed B05441 wrist (45mm) is not compatible with **size 6** hands. Recommended use is with friction housing (p.143) and hand plate (p.143).

Part Number	Size	UK Size	Left/Right	Application	Wrist
SP592L	5 ½	2" x 3"	Left	Child 2-5 years	40mm (1 ½")
SP592R	5 ½	2" x 3"	Right	Child 2-5 years	40mm (1 ½")
SP589L	6	2 ¼" x 4"	Left	Child 5-8 years	45mm (1 ¾")
SP589R	6	2 ¼" x 4"	Right	Child 5-8 years	45mm (1 ¾")



### Spare Parts Operating Cable

Spare front pull operating cable cord for Child Cable Operated Hands.

Part Number	Hand Size	UK Size
B07365	5 ½, 6	2", 2 ¼"