

Wrist Units - Child Wrist, Fixed Wrist and Other Hand Plates



Child 2-Knob Wrist Unit

Part Number	Wrist Size	Description	Weight	Height
1 B05143	40mm (1 1/2")	Quick release, lock, 12 position rotation	57g (2 oz)	31mm (1 1/4")
B05142	45mm (1 3/4")	Quick release, lock, 12 position rotation	84g (3 oz)	31mm (1 1/4")
NSC0375	40mm (1 1/2") + 45mm (1 3/4")	Replacement Screws		

Supplied with parallel housing and three screws.

Child 2-Knob Wrist Hand Plate

Part Number	Wrist Size	Description	Weight
2 B07030	40mm (1 1/2" x 1/2" x 3/8")	1 1/2" Child 2-Knob Rotary	26g (1 oz)
B07029	45mm (1 3/4" x 1/2" x 3/8")	1 3/4" Child 2-Knob Rotary	26g (1 oz)

Child 2-Knob Wrist Adapter Plate

Part Number	Wrist Size	Internal Thread	Weight
3 B00897	40mm (1 1/2" x 1/2" x 3/8")	5/16" BSF (UK Juvenile Standard)	22g (1oz)
B00821	45mm (1 3/4" x 1/2" x 3/8")	5/16" BSF (UK Juvenile Standard)	27g (1oz)
B00905	45mm (1 3/4" x 1/2" x 3/8")	3/8" BSW (UK Adult Standard)	45g (1 1/2 oz)

Fixed Wrist Connector and Former - Round

Part Number	Wrist Size	Former
4 B05963	40mm (1 1/2")	4a SB12136/1
B05441	45mm (1 3/4")	SB12136/2
B05440	50mm (2")	SB12136/3
B05439	54mm (2 1/8")	SB12136/4

Fixed Wrist Connector and Former - Oval

Part Number	Wrist Size	Former
5 B06252	45mm (1 3/4")	5a B06986
B06256	50mm (2")	B06985
B06255	54mm (2 1/8")	B06984

Hand Plate - North American - 1/2"

Part Number	Wrist Size	External Thread	Shape	Weight
6 B10530	40mm (1 1/2" x 1/16" x 1/2" 20 TPI)	1/2" x 20 TPI	Round - Steeper Hand Attachment	10g (1/2 oz)
B22758	40mm (1 1/2" x 1/16" x 1/2" 20 TPI)	1/2" x 20 TPI	Round - Otto Bock Hand Attachment	16g (1/2 oz)
B22757	45mm (1 3/4" x 1/16" x 1/2" 20 TPI)	1/2" x 20 TPI	Round - Steeper Hand Attachment	16g (1/2 oz)
B22756	50mm (2" x 1/16" x 1/2" 20 TPI)	1/2" x 20 TPI	Round - Steeper Hand Attachment	17g (1/2 oz)
B22755	54mm (2 1/8" x 1/16" x 1/2" 20 TPI)	1/2" x 20 TPI	Round - Steeper Hand Attachment	19g (3/5 oz)

Both B22757, B22756 and B22755 are also compatible with Ottobock fastening arrangement.

7 B06257	45mm (1 3/4" x 1/16" x 1/2" 20 TPI)	1/2" x 20 TPI	Oval	10g (3/5 oz)
----------	-------------------------------------	---------------	------	--------------

Note: For use with Foam Hand with Oval base (2 1/2" - 3 3/4").

Hand Plate - European - M12

Part Number	Wrist Size	External Thread	Shape	Weight
8 B22768	40mm (1 1/2" x 3/4" x M12)	M12 x 1.5	Round	15g (1/2 oz)
B07786	45mm (1 3/4" x 3/4" x M12)	M12 x 1.5	Round	15g (1/2 oz)
B06942	50mm (2" x 3/4" x M12)	M12 x 1.5	Round	15g (1/2 oz)

Both B06942 and B07786 are also compatible with Ottobock fastening arrangement.