



Technical Information: Cable Operated and Spring Hands

Size	Width at Knuckles	Cable Pull	Max Opening	Weight	Base to Thumb Tip
5 ½	50mm (2")	29mm (1 ⅛")	44mm (1 ¾")	84g (3oz)	76mm (3")
6	57mm (2 ¼")	25mm (1")	50mm (2")	130g (4 ½oz)	102mm (4")
6 ¾	63mm (2 ½")	60mm (2 ⅜")	58mm (2 ¼")	224g (8oz)	108 (4 ¼")
7 ¼	70mm (2 ¾")	40mm (1 ⅝")	90mm (3 ½")	323g (11oz)	133mm (5 ¼")
7 ¾	82mm (3 ¼")	40mm (1 ⅝")	90mm (3 ½")	342g (12oz)	133mm (5 ¼")
8 ¼	87mm (3 ½")	42mm (1 ⅝")	95mm (3 ¾")	366g (13oz)	137mm (5 ⅜")

Child Cable Operated Hand Front Pull - Voluntary Opening

Part Number	Size	UK Size	Left / Right	Application	Wrist	Glove No.
1 B08361	5 ½	2" x 3"	Left	Child 2-5 years	40mm (1 ½")	46
B08360	5 ½	2" x 3"	Right	Child 2-5 years	40mm (1 ½")	47
C07331	6	2 ¼" x 4"	Left	Child 5-8 years	45mm (1 ¾")	44
C07330	6	2 ¼" x 4"	Right	Child 5-8 years	45mm (1 ¾")	45

Child Cable Operated Hand Back Pull - Voluntary Opening

Part Number	Size	UK Size	Left / Right	Application	Wrist	Glove No.
2 B08363	5 ½	2" x 3"	Left	Child 2-5 years	40mm (1 ½")	46
B08362	5 ½	2" x 3"	Right	Child 2-5 years	40mm (1 ½")	47
C07931	6	2 ¼" x 4"	Left	Child 5-8 years	45mm (1 ¾")	44
C07932	6	2 ¼" x 4"	Right	Child 5-8 years	45mm (1 ¾")	45

Child Spring Hand

Part Number	Size	UK Size	Left / Right	Application	Wrist	Glove No.
3 SP592L	5 ½	2" x 3"	Left	Child 2-5 years	40mm (1 ½")	46
SP592R	5 ½	2" x 3"	Right	Child 2-5 years	40mm (1 ½")	47
SP589L	6	2 ¼" x 4"	Left	Child 5-8 years	45mm (1 ¾")	44
SP589R	6	2 ¼" x 4"	Right	Child 5-8 years	45mm (1 ¾")	45

Child Cable Operated Hand - Spare Parts Operating Cable

Part Number	Description	Size	UK Size
4 B07365	Operating Cable, Front Pull	5 ½, 6	2", 2 ¼"
B10734	Operating Cable, Back Pull	5 ½, 6	2", 2 ¼"