


Electrode important information


Keep this card for reference




1. Do not immerse the electrode in water.



2. Do not expose the electrode to dirt, dust or oil.



3. Do not attempt to disassemble the electrode.



4. Do not expose the electrode to shock or vibration.

Electrode important information

Keep this card for reference

- 1. Ensure the electrode contacts are regularly cleaned. Use a lightly moistened cloth with soapy water.**
- 2. Never clean the electrode with a solvent based cleanser.**
- 3. If your prosthesis stops working contact your prosthetist.**
- 4. If your prosthesis becomes uncomfortable or the control becomes unpredictable contact your prosthetist.**

Steeper Group,
Unit 3 Stourton Link,
Intermezzo Drive,
Leeds,
LS10 1DF

RSLIT328 Issue 3 © 2017 Steeper

