



Heavy Duty Wrist Unit				
Part Number	Size	Description	Weight	Height
1 SP249	50mm (2")	Quick Release, lock, 12 position rotation	70g (2 ½ oz)	27mm (1")
B01155	54mm (2 ⅛")	Quick Release, lock, 12 position rotation	72g (2 ½ oz)	27mm (1")
SP136	58mm (2 ½")	Quick Release, lock, 12 position rotation	74g (2 ½ oz)	27mm (1")

Wrist Housing - Parallel				
Part Number	Size	Description	Weight	Height
2 B00573	50mm (2")	Parallel, bonded into forearm	20g (½ oz)	25mm (1")
B00572	54mm (2 ⅛")	Parallel, bonded into forearm	20g (½ oz)	25mm (1")





Wrist Housing - Tapered				
Part Number	Size	Description	Weight	Height
3 B00825	50mm (2")	3° Taper, Laminated into forearm	27g (1oz)	38mm (1 ½")
B00827	54mm (2 ⅛")	3° Taper, Laminated into forearm	28g (1oz)	38mm (1 ½")
B00826	54mm (2 ⅛")	6° Taper, Laminated into forearm	28g (1oz)	38mm (1 ½")
B00828	58mm (2 ½")	5° Taper, Laminated into forearm	28g (1oz)	38mm (1 ½")
B01017	58mm (2 ½")	7° Taper, Laminated into forearm	28g (1oz)	38mm (1 ½")

Hand Connection Plate				
Part Number	Size	Description	Weight	Height
4 B01179	50mm (2" x ¾" x ½")	2" Heavy Duty Wrist	56g (2 oz)	31mm (1 ¼")
B01172	54mm (2 ⅛" x ¾" x ½")	2 ⅛" and 2 ½" Heavy Duty Wrist	56g (2 oz)	31mm (1 ¼")

**Note:** Both B01172 and B01179 are now also compatible with Ottobock fastening arrangement.

Adapter Plates - Flat					
Part Number	Size	Description	Compatible With	Weight	Height
5 B00967	50mm (2" x ¾" x ½")	⅜" BSW (UK Thread Internal)	2" Heavy Duty Wrist	46g (1 ½ oz)	34mm (1 ⅜")
B00900	54mm (2 ⅛" x ¾" x ½")	⅜" BSW (UK Thread Internal)	2 ⅛" , 2 ½" Heavy Duty Wrist	46g (1 ½ oz)	34mm (1 ⅜")
6 B06568	50mm (2" x ¾" x ½")	½" x 20TPI (US Thread Internal)	2" , 2 ⅛" , 2 ½" Heavy Duty Wrist	46g (1 ½ oz)	49mm (2")
B22264	54mm (2 ⅛" x ¾" x ½")	M12 (European Thread Internal)	2 ⅛" , 2 ½" Heavy Duty Wrist	46g (1 ½ oz)	54mm (2 ⅛")

Adapter Plates - Domed					
Part Number	Size	Description	Compatible With	Weight	Height
7 B00965	54mm (2 ⅛" x ¾" x ½")	⅜" BSW (UK Thread Internal) Allows free rotation with lock feature	2 ⅛" , 2 ½" Heavy Duty Wrist	55g (2 oz)	34mm (1 ⅜")
7a B00966	54mm (2 ⅛" x ¾" x ½")	⅜" BSW (UK Thread Internal) Allows free rotation without lock feature	2 ⅛" , 2 ½" Heavy Duty Wrist	55g (2 oz)	34mm (1 ⅜")



Wrist Units - Heavy Duty Wrist Spare Parts

Heavy Duty Wrist Unit - Spare Parts

Part Number	Size	Description	Supplied
1 E17985	50mm, 54mm (2" , 2 1/8" ) Heavy Duty Wrists	Operating Knob Kit 2" , 2 1/8"	Each
E17987	58mm (2 1/2") Heavy Duty Wrists	Operating Knob Kit 2 1/2"	Each
2 B05435	All Heavy Duty Wrists	Sliding Cover - Parallel	Each
B05368	All Heavy Duty Wrists	Sliding Cover - Tapered	Each
3 NPN0217	All Heavy Duty Wrists	Locking Mills Pin	Each
4 B00748	All Heavy Duty Wrists	Rubber Locking Ring	Each
5 B00744	All Heavy Duty Wrists	Snap Catch Blade	Each
6 B00746	All Heavy Duty Wrists	Catch Blade Spring	Each
7 NPN0007	All Heavy Duty Wrists	Split Pin	Each
8 B00745	All Heavy Duty Wrists	Plug	Each

Heavy Duty Wrist Unit - Tools

Part Number	Size	Description	Supplied
9 B18631	All Heavy Duty Wrists	Assembly Removal Tool	Each