




Friction Wrist				
Part Number	Wrist Size	Application	Weight	Height
1 C19011	40mm (1 1/2")	Child	31g (1oz)	31mm (1 1/4")
C12848	45mm (1 3/4")	Child / Juvenile	37g (1oz)	31mm (1 1/4")
C12810	50mm (2")	Adolescent / Adult	49g (1 3/4 oz)	31mm (1 1/4")
C12824	54mm (2 1/8")	Adult	72g (2 1/2 oz)	31mm (1 1/4")

Friction Wrist Spare Parts				
Part Number	Size	Description	Supplied	
2 C18132	40mm (1 1/2")	Friction Ring	Each	
B12852	45mm (1 3/4")	Friction Ring	Each	
B12815	50mm (2")	Friction Ring	Each	
B12821	54mm (2 1/8")	Friction Ring	Each	
2a B18133	40mm (1 1/2")	Protective Rings	Each	
B12850	45mm (1 3/4")	Protective Rings	Each	
B12813	50mm (2")	Protective Rings	Each	
B12819	54mm (2 1/8")	Protective Rings	Each	

Friction Wrist Adapter Plate				
Part Number	Wrist Size	Internal Thread	Weight	Height
3 B18161	40mm (1 1/2")	5/16" BSF (UK Juvenile Standard)	12g (1/2 oz)	14.2mm (1/2")
B22959	40mm (1 1/2")	M12 x 1.5 (European Standard)	11g (1/2 oz)	18mm (3/4")
B19639	40mm (1 1/2")	1/2" x 20 TPI (US Standard)	10g (1/2 oz)	18mm (3/4")
B18169	45mm (1 3/4")	5/16" BSF (UK Juvenile Standard)	13g (1/2 oz)	18mm (3/4")
B18170	45mm (1 3/4")	3/8" BSW (UK Adult Standard)	13g (1/2 oz)	18mm (3/4")
B18171	45mm (1 3/4")	M12 x 1.5 (European Standard)	13g (1/2 oz)	18mm (3/4")
B18172	45mm (1 3/4")	1/2" x 20 TPI (US Standard)	13g (1/2 oz)	18mm (3/4")
B18173	50mm (2")	3/8" BSW (UK Adult Standard)	15g (1/2 oz)	18mm (3/4")
B18174	50mm (2")	M12 x 1.5 (European Standard)	15g (1/2 oz)	18mm (3/4")
B18175	50mm (2")	1/2" x 20 TPI (US Standard)	15g (1/2 oz)	18mm (3/4")
B18176	54mm (2 1/8")	3/8" BSW (UK Adult Standard)	19g (1/2 oz)	18mm (3/4")
B18177	54mm (2 1/8")	M12 x 1.5 (European Standard)	19g (1/2 oz)	18mm (3/4")
B18178	54mm (2 1/8")	1/2" x 20 TPI (US Standard)	19g (1/2 oz)	18mm (3/4")

Friction Housing				
Part Number	Description	Thread	Weight	Height
4 NHD55108	Friction Housing 40mm (1 1/2")	Internal Thread 1/2" UNF x 20 TPI	23g (1oz)	16mm (5/8")
NHD55109	Friction Housing 45mm (1 3/4")	Internal Thread 1/2" UNF x 20 TPI	28g (1oz)	16mm (5/8")
NHD55110	Friction Housing 50mm (2")	Internal Thread 1/2" UNF x 20 TPI	35g (1oz)	16mm (5/8")
NHD55108-M12	Friction Housing 40mm (1 1/2")	Internal Thread M12	23g (1oz)	16mm (5/8")
NHD55109-M12	Friction Housing 45mm (1 3/4")	Internal Thread M12	28g (1oz)	16mm (5/8")
NHD55110-M12	Friction Housing 50mm (2")	Internal Thread M12	35g (1oz)	16mm (5/8")



Friction Wrist Hand Connection Plates (No Inserts)

Part Number	Wrist Size	Weight	Height
5 B18165	40mm (1 1/2")	4g (1/8 oz)	8.4mm (3/8")
B18215	45mm (1 3/4")	5g (1/8 oz)	8.9mm (3/8")
B12811	50mm (2")	14g (1/2 oz)	8.9mm (3/8")
B12817	54mm (2 1/8")	17g (5/8 oz)	8.9mm (3/8")

Friction Wrist Hand Connection for Myoelectric, Spring and Cable Operated Hands

Part Number	Wrist Size	Description	Weight
6 KFW045	45mm (1 3/4")	Friction Wrist handplate kit for above hands	12g (1/2 oz)
KFW050	50mm (2")	Friction Wrist handplate kit for above hands	14g (1/2 oz)
KFW054	54mm (2 1/8")	Friction Wrist handplate kit for above hands	17g (1/2 oz)

Short Wrist

Part Number	Size	Description	Weight	Height
7 BBSW50	50mm (2")	1/2" x 20 TPI Internal Thread	22g (3/4 oz)	10mm (3/8")
BBSW45	45mm (2")	1/2" x 20 TPI Internal Thread	22g (3/4 oz)	10mm (3/8")
BBSW50-M12	50mm (2")	M12 Internal Thread	22g (3/4 oz)	10mm (3/8")
7a B27336	50mm (2")	1/2" x 20 TPI Internal Thread - No Lamination Collar	20g (3/4 oz)	10mm (3/8")

Wooden Wrist

Part Number	Size	Description	Weight	Height
8 N621020100	60mm (2 3/8")	M12 thread. Parallel, bonded into forearm	42.5g (1 1/2")	20mm (3/4")