

Modular Arm System - System Components

Upper Limb Prosthetics
Ordering Information
Myoelectric Hands
Myoelectric Elbows
Myoelectric Power & Inputs
Cable Operated & Spring Hands
Functional Devices
Cosmetic Hands & Gloves
Modular Arm Systems
Wrist Units
Mech. Elbow Units & Joints
Suspension & Control Systems
Materials & Tools



Friction Wrist Adapter

For use with 50mm (2") or 54mm (2 1/8") Friction Wrist and Connection Plates, see pages 141 - 143 for more details. Supplied with Split Pin NPN009.


Part Number	Description	Weight	Height
B17951	Multi-Axial Wrist Housing 50mm (2")	41g (1 1/2 oz)	32mm (1 1/4")
B17950	Multi-Axial Wrist Housing 54mm (2 1/8")	37g (1 1/2 oz)	32mm (1 1/4")



Multi-Axial Wrist Housing

Multi-Axial Wrist Housing for attachment to a 54mm Hand Connection Plate. Ball Insert supplied separately (see below).

Part Number	Description	Weight	Height
B22505	Multi-Axial Wrist Housing 54mm (2 1/8")	51g (1 1/2 oz)	25mm (1")



Flexion Ball Inserts

Part Number	Description	Thread	Weight	Height
1 B22763	Free Flexion Ball Insert - Allows rotation, flexion and extension	Internal Thread 1/2" UNF x 20 TPI	18g (1/2 oz)	24mm (1")
B22503	Free Flexion Ball Insert - Allows rotation, flexion and extension	Internal Thread 3/8" BSW	20g (1/2 oz)	24mm (1")
2 B22764	Locked Flexion Ball Insert - Allows rotation only	Internal Thread 1/2" UNF x 20 TPI	20g (1/2 oz)	27mm (1")
B22765	Locked Flexion Ball Insert - Allows rotation only	Internal Thread 3/8" BSW	21g (1/2 oz)	26mm (1")