



## Quick Disconnect Wrist Unit

Part Number	Description	Hand Size	Height	Weight
5 QDAWA-40	40mm Wrist Assembly. Includes Lamination Ring, Coupling Piece and Coaxial Plug	6 ¾	20mm (¾")	38g (1 ⅜ oz)
QDAWA-45	45mm Wrist Assembly. Includes Lamination Ring, Coupling Piece and Coaxial Plug	7 ¼	25mm (1")	38g (1 ⅜ oz)
QDAWA-50	50mm Wrist Assembly. Includes Lamination Ring, Coupling Piece and Coaxial Plug	7 ¾	25mm (1")	42g (1 ½ oz)
QDAWA-54	54mm Wrist Assembly. Includes Lamination Ring, Coupling Piece and Coaxial Plug	8 ¼	25mm (1")	44g (1 ½ oz)

## Quick Disconnect Wrist - Spares

Part Number	Description	Hand Size	Height
6 QDALR-40	40mm Lamination Ring	6 ¾	20mm (¾")
QDALR-45	45mm Lamination Ring	7 ¼	23mm (⅞")
QDALR-50	50mm Lamination Ring	7 ¾	23mm (⅞")
QDALR-54	54mm Lamination Ring	8 ¼	23mm (⅞")
7 QDACP	Coupling Piece. Suitable for 40, 45, 50 and 54mm Quick Disconnect Wrist. Includes lock ring.		
8 QDACO	Coaxial Plug. Suitable for 40, 45, 50 and 54mm Quick Disconnect Wrist. Includes lock ring and screw.		

## Battery

Part Number	Description	Chemistry	Capacity	Dimensions	Weight
9 BLPA72	7.2V Lithium Polymer Battery	Lithium Polymer	800mAh	68 x 33 x 17 mm (2 ⅝" x 1 ¼" x ⅝")	45g (1 ½ oz)

## Battery Charger

Part Number	Description	Region
10 LPDC01	Lithium Polymer Battery Desktop Charger	World

## Battery Receptacles

Forearm mounted battery receptacles are supplied in kits with all necessary formers

Part Number	Description	Battery
11 BHQDA72	Lithium Polymer Battery Holder with 2-Pin Output	BLPA72

**Note:** Suitable for connection to Quick Disconnect Wrist connection

12 BHA72	Lithium Polymer Battery Holder with Locking Connector	BLPA72
----------	---	--------

**Note:** Suitable for connection to Connector Block Cable CBBHA72.