



Durable, yet lightweight.

Designed with the basic principles of durability and ease of use firmly in mind.

Using high strength aerospace alloy inside a flexible PVC shell, these prosthetic hands are also compatible with our range of full-length Elegance™ cosmetic gloves in 19 lifelike TrueFinish™ shades for enhanced realism.

Our advanced materials and manufacturing processes ensure that these cable and spring operated prosthetic hands are detailed, durable and lightweight options for patients.

Featuring:

- Body-powered cable control or manual operation
- Voluntary opening - front pull, back pull or spring close options
- 3 cable terminations supplied
- Four sizes available from child to adult
- Plain hand base, M12 x 1.5 or ½ x 20 TPI threaded stem



Cable and Spring Operated Mechanical Hands



To Order

Prefix the number shown in the table with TG for TrueFinish PVC or TSG for TrueFinish Silicone.

Then, add the preferred skin shade. E.g. **TSG300/E4** orders a 8¼ left adult silicone glove in TrueFinish colour E4.

To order a color swatch use part number **E24655**.



Ordering Information - Spring Operated Hands								
Part No. Plain Base	Part No. M12 x 1.5 Stud	Part No. ½ x 20 TPI Stud	Size	Left / Right	Application	Wrist Diameter	TG = TrueFinish PVC TSG = TrueFinish Silicone	
							Male Glove	Female Glove
SPA675L	SPA675LEU	SPA675LNA	6 ¾	Left	Child (Approx.8-12)	45mm	310	310
SPA675R	SPA675REU	SPA675RNA	6 ¾	Right	Child (Approx.8-12)	45mm	311	311
SPA725L	SPA725LEU	SPA725LNA	7 ¼	Left	Adolescent / Adult	50mm	306	304
SPA725R	SPA725REU	SPA725RNA	7 ¼	Right	Adolescent / Adult	50mm	307	305
SPA775L	SPA775LEU	SPA775LNA	7 ¾	Left	Adult	54mm	302	312
SPA775R	SPA775REU	SPA775RNA	7 ¾	Right	Adult	54mm	303	313
SPA825L	SPA825LEU	SPA825LNA	8 ¼	Left	Adult	54mm	300	
SPA825R	SPA825REU	SPA825RNA	8 ¼	Right	Adult	54mm	301	

Ordering Information - Cable Operated Hands - Front Pull Voluntary Opening								
Part No. Plain Base	Part No. M12 x 1.5 Stud	Part No. ½ x 20 TPI Stud	Size	Left / Right	Application	Wrist Diameter	TG = TrueFinish PVC TSG = TrueFinish Silicone	
							Male Glove	Female Glove
MHA675LF	MHA675LFEU	MHA675LFNA	6 ¾	Left	Child (Approx.8-12)	45mm	310	310
MHA675RF	MHA675RFEU	MHA675RFNA	6 ¾	Right	Child (Approx.8-12)	45mm	311	311
MHA725LF	MHA725LFEU	MHA725LFNA	7 ¼	Left	Adolescent / Adult	50mm	306	304
MHA725RF	MHA725RFEU	MHA725RFNA	7 ¼	Right	Adolescent / Adult	50mm	307	305
MHA775LF	MHA775LFEU	MHA775LFNA	7 ¾	Left	Adult	54mm	302	312
MHA775RF	MHA775RFEU	MHA775RFNA	7 ¾	Right	Adult	54mm	303	313
MHA825LF	MHA825LFEU	MHA825LFNA	8 ¼	Left	Adult	54mm	300	
MHA825RF	MHA825RFEU	MHA825RFNA	8 ¼	Right	Adult	54mm	301	

Ordering Information - Cable Operated Hands - Back Pull Voluntary Opening								
Part No. Plain Base	Part No. M12 x 1.5 Stud	Part No. ½ x 20 TPI Stud	Size	Left / Right	Application	Wrist Diameter	TG = TrueFinish PVC TSG = TrueFinish Silicone	
							Male Glove	Female Glove
MHA675LB	MHA675LBEU	MHA675LBNA	6 ¾	Left	Child (Approx.8-12)	45mm	310	310
MHA675RB	MHA675RBEU	MHA675RBNA	6 ¾	Right	Child (Approx.8-12)	45mm	311	311
MHA725LB	MHA725LBEU	MHA725LBNA	7 ¼	Left	Adolescent / Adult	50mm	306	304
MHA725RB	MHA725RBEU	MHA725RBNA	7 ¼	Right	Adolescent / Adult	50mm	307	305
MHA775LB	MHA775LBEU	MHA775LBNA	7 ¾	Left	Adult	54mm	302	312
MHA775RB	MHA775RBEU	MHA775RBNA	7 ¾	Right	Adult	54mm	303	313
MHA825LB	MHA825LBEU	MHA825LBNA	8 ¼	Left	Adult	54mm	300	
MHA825RB	MHA825RBEU	MHA825RBNA	8 ¼	Right	Adult	54mm	301	

SteeperUSA, 8666 Huebner Road, Suite 112, San Antonio, TX 78240, USA

www.steeperusa.com

STPPRO38 Issue 2 2020

